


## STL Export Instructions (myitero.com)

1. Open Google Chrome or Microsoft Edge (recommended browsers).




2. Go to [www.MyiTero.com](http://www.MyiTero.com), you must login with the same email and password that you use on the scanner.


3. Upon logging in, go to "Patients."



4. Once you have located your patient's scan, click "Export".




5. A window will pop-up for you to choose the CAD/CAM program the lab uses. Mabel recommends 3Shape or if the program does not appear, select the “Implants and other exports-Open Model” option.



6. Once you choose the CAD/CAM program, you can click on Export.


6.1 You will see a download icon on the upper right side of the screen.



6.2 Wait until the download is completed.

7. Once the folder is downloaded, you can open the files in the folder where they have been saved by clicking the arrow as shown below.

7.1 Chrome (browser)



## 7.2 Edge (browser)



8. Once you show in folder, you will see a zipped folder under the name of "OrthoCAD\_Export\_99999999 (the order ID of the scan will be shown)". You can send that compressed folder to the lab or the person who will work on the scans via email.

OrthoCAD\_Export\_99999999      10/23/2018 10:23 ...      Compressed (zipp...      16,796 KB

Please do not hesitate to contact us at 877.622.3533 or email us at [support@mabeldental.com](mailto:support@mabeldental.com) or contact iTero directly at 800-577-8767, or via email: [iterosupport@aligntech.com](mailto:iterosupport@aligntech.com) if you have any other concerns.

**If the file is too large to email, please see the next page for file upload instructions:**

## Are your digital files too large to email?

Please transfer your large digital scan files to us by following these instructions:


We recommend the safe and secure website WeTransfer.com and here's how to use it:

1. Go to [WeTransfer.com](https://www.wetransfer.com).
2. Choose the "Take Me to Free" option.
3. Choose "I Agree" to the Terms and Condition
4. Click "+ Add your files" and select the files from your computer that you want to transfer. We recommend to also include the Prescription Rx.
5. In the "Email to" field, enter support@mabeldental.com
6. Enter your email address in the "Your email" field.
7. Enter an optional message.
8. Click the "Transfer" button.  
Voila! The file(s) has been sent!

If you did not include the Prescription Rx, we recommend that you fax it to 234.678.5341 or email it to: support@mabeldental.com.

You may view all of our Rx [here](#) on our website.

You'll receive confirmation from us within 1 business day.



The screenshot shows the WeTransfer.com interface. At the top, there is a blue circular button with a white plus sign and the text "Add your files". Below this, there are three input fields: "Email to", "Your email", and "Message". At the bottom, there is a grey button with the text "Transfer" and a small circular icon with three dots to its left.

\*please note this may not reflect the exact look of the screen



877.622.3533 | mabeldental.com