

# Digital Denture Printed Try-In Checklist

## Ensure a Perfect Fit for Your Patient

The Digital Denture Printed Try-In Patient Checklist helps you and the patient evaluate the fit, function, and esthetics of the digitally printed monochrome try-in. Include the completed form to keep your case moving smoothly and to ensure precision in the final printed denture.

Doctor Name:

Practice / Lab Name:

Patient Name:

Date:

Return Case by 5pm:

## Overview

- The printed try-in is an exact replica of how the final denture will be including fit, midline and teeth size.
- Adjust & equilibrate the digital printed try-in. Then scan the adjusted try-in denture and send it to our lab.
- Provide exacting adjustment instructions expressed in millimeters of change.
- In addition to the checklist comments, if desired, express adjustment instructions with a bur on the try-in.

## Step-By-Step Checklist

**1. Insert the try-in in the patient's mouth and adjust sharp spots or pressure areas.**

**2. Does the printed try-in have good retention?** ☐ Yes ☐ No

If no, please take a new impression of the try-in.

Comments: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

**3. Is the Vertical Dimension of Occlusion (VDO) correct?** ☐ Yes ☐ No

If no, adjust the occlusion and make a new record.

Comments: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

## Step-By-Step Checklist *continued*

### 4. Are there areas of over-extension? ☐ Yes ☐ No

If yes, grind away any over-extension and take a new impression to send to our lab.

Comments: \_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

### 5. Is the midline in position? ☐ Yes ☐ No

If no, mark adjusted midline on try-in denture with thin diameter bur and take a new impression.

Comments: \_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

### 6. Is the lip support adequate? ☐ Yes ☐ No

If no, indicate facial or lingual – indicate adjustment +/- in mm. Evaluate the upper teeth on the wet/dry line of the lower lip. Evaluate “SH” sound for lower teeth.

Comments: \_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

### 7. Is the incisal edge position correct? ☐ Yes ☐ No

If no, determine if teeth are too high or low – indicate adjustment +/- in mm. Check “F” and “V” sounds.

Comments: \_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

### 8. Is there a cant? ☐ Yes ☐ No

If yes, provide teeth adjustment, left or right – indicate adjustment +/- in mm.

Comments: \_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

### 9. Check the cervical of #6-#11 for proper placement. ☐ Yes ☐ No

If incorrect, indicate in millimeters how much more/less tissue/tooth needs to show.

Comments: \_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

### 10. Select final printed denture tooth shade:

☐ A1 ☐ A2 ☐ A3 ☐ B1 ☐ Bleach

Submit this checklist for noted  
adjustments and final design.

### 11. How would you like to proceed?

☐ Return for another printed try-in ☐ Return a final printed denture