

Your Smile Partner



MARCH 2025 LAB TO LAB PRICE LIST

877.622.3533 | mabeldental.com





PORCELAIN FUSED TO METAL

| 8 in-lab days 6 in-lab days (digital) | |
|---|------------|
| Non-precious | 49 |
| Non-precious coping only | 25 |
| Noble | 35 + alloy |
| with metal occlusal | 60 + alloy |
| High Noble (White) | 35 + alloy |
| with metal occlusal | 74 + alloy |
| High Noble (Yellow) | 35 + alloy |
| with metal occlusal | 87 + alloy |

CUSTOM ABUTMENTS

| 8 in-lab days 6 in-lab days (digital) | |
|---|-----|
| Custom Abutment – Cast* 50 + gold | |
| Custom Abutment - Titanium | 119 |
| Custom Abutment - Zirconia | 139 |

Price above includes abutment only. Additional implant parts are not included.*Requires UCLA abutment, not included.

RESTORATIONS OVER IMPLANTS

| 8 in-lab days 6 in-lab days (digital) | |
|---|------------|
| Full Contour Zirconia | 79 |
| IPS e.max® layered | 89 |
| Multi-Layered Zirconia | 85 |
| Porcelain Fused to Zirconia | 85 |
| PFM – Non-precious crown | 69 |
| PFM – Semi-precious crown | 55 + alloy |
| PFM – High Noble crown | 55 + alloy |
| Full-Cast | add \$20 |

Implant crown prices above include only model work and final crown for cemented and screw-retained crowns.

Additional labor and implant parts are not included.

| 25 |
|----|
| 25 |
| 10 |
| 30 |
| |

DIAGNOSTIC WAX-UP

| 8 in-lab days 6 in-lab days (digital) | |
|---|----|
| Diagnostic Wax-Up (per unit) | 20 |

ZIRCONIA

| 8 in-lab days 6 in-lab days (digital) | | |
|---|----|--|
| Full Contour coping only 59 39 | | |
| Multi-Layered | 65 | |
| Porcelain Fused to Zirconia | 65 | |
| Temporary PMMA (monolithic) | 25 | |

IPS E.MAX®

| 8 in-lab days 6 in-lab days (digital) | |
|---|----|
| IPS e.max® layered | 69 |
| IPS e.max Inlay/Onlay | 69 |
| IPS e.max Veneer | 69 |

FULL-CAST GOLD

| 8 in-lab days 6 in-lab days (digital) | |
|---|------------|
| Non-precious | 49 |
| 2% Yellow Gold | 30 + alloy |
| 40% Yellow Gold | 30 + alloy |
| 40% White Gold | 30 + alloy |
| 60% Yellow Gold | 30 + alloy |

MISCELLANEOUS

| All Porcelain Butt Shoulder 360° | 10 |
|--------------------------------------|--------|
| Butt Shoulder 180° (Porcelain) | 10 |
| Cingulum Guideplane / Rest | 7 |
| Contour Crown | 10 |
| Crown to Existing Partial (per unit) | 30 |
| Crown to Future Partial (per unit) | 40 |
| Digital Design only (per unit) | 10 |
| Guideplane / Rest Seat | 7 |
| Printed Model – Quad Full (arch) | 8 12 |
| Post & Core - Non-precious | 20 |
| Post & Core - Noble | 55 |
| Shade Change Fee | 25 |
| Transfer / Reduction Coping | 14 |

Turnaround Time

Physical Cases: 8 days | Digital Cases: 6 days add 2 days for implants

Digital Design: 24 - 48 hours

| FULL DENTURES | | |
|--|-------------|--|
| STANDARD COMPLETE | 99 | |
| includes mounting, set-up, teeth, festooning/stippling, Lucitone 199® finish & patient ID | | |
| PREMIUM COMPLETE | 157 | |
| includes mounting, set-up, premium teeth, festooning/stippling, Lucitone 199 finish & patient ID | | |
| Set-up | 49 | |
| Finish | 50 | |
| Premium Teeth – Anterior | 5 per tooth | |
| Premium Teeth – Posterior | 3.5 / tooth | |

| METAL PARTIALS | |
|---|-----|
| SLM COMPLETE | 149 |
| includes SLM frame, mounting, set-up, teeth, Lucitone 199 finish & patient ID | |
| VITALLIUM 2000 COMPLETE | 169 |
| Includes Vitallium 2000 frame, mounting, setup, teeth, Lucitone 199 finish & patient ID | |
| Add Flexible Clasp or Tooth Colored Clasp | 25 |
| Frame only – SLM (selective laser melting) | 69 |
| Frame only – Titanium | 79 |
| Frame only - Vitallium® 2000 (cast) | 89 |
| Set-up – Metal Partial | 40 |
| Finish – Metal Partial Lucitone 199 | 40 |

| FLEXIBLE PARTIALS | | |
|---|-----|--|
| COMPLETE includes mounting, set-up, teeth & TCS [®] finish | 129 | |
| Unilateral/Nesbit (3 teeth max) | 99 | |
| Set-up – Flexible Partial | 50 | |
| Finish – Flexible Partial (TCS) | 79 | |
| Add Flexible saddles/clasps to frame (per unit) | 25 | |

| ACRYLIC FLIPPERS | |
|--|-----|
| Acrylic Flipper 1 - 4 teeth | 79 |
| Acrylic Flipper 5-10 teeth (excludes teeth) | 99 |
| Acrylic Flipper 11+ teeth (excludes teeth) | 109 |
| Add wrought wire clasps | 20 |
| includes mounting, set-up, teeth, finish & ID per arch | |
| | |

VISICLEAR™ PARTIALS COMPLETE 169 includes VisiClear™ frame, mounting, set-up, teeth, Lucitone 199 finish & patient ID Frame only - VisiClear 89 Set-up – VisiClear Partial 40 Finish – VisiClear Partial Lucitone 199 40

| NIGHTGUARDS / BITE SPLINTS | |
|---------------------------------|----|
| Hard / Soft Combo (vacuum form) | 65 |
| Hard / Soft – Milled Impak™ | 65 |
| Hard - Processed Acrylic | 59 |
| Hard – 3D Printed | 49 |
| Soft (vacuum form) | 35 |

| DURACETAL® PARTIALS | | |
|---|-----|--|
| COMPLETE Includes Duracetal frame, dup model, mounting, set-up, teeth, Lucitone 199 finish & patient ID | 179 | |
| Frame only – Duracetal | 99 | |
| Set-up – Duracetal Partial | 40 | |
| Finish – Duracetal Partial | 40 | |
| Available Shades - A1, A2, A3.5 & B1 | | |

DIGITAL DENTURES - 3D Printed

COMPLETE
includes design, set-up, teeth, finish & ID
must be straight to finish, no try-ins

SNAP-ON SMILE®

| Snap-On Smile we are an authorized laboratory | 299 |
|---|-----|
| Available Shades - A1, A2, A3.5, B1 & Blea | ach |
| Smile Temps – vacuum form, economy version | 149 |
| Available Shades - A1, A2, OM1 | |

| REPAIRS | |
|---------------------------------|----|
| Acrylic Repair (1 tooth) | 35 |
| each additional tooth | 7 |
| Fracture Repair | 40 |
| Rebase / Jump | 50 |
| Reline - Hard | 50 |
| Reline with soft liner | 80 |
| Laser Weld | 30 |
| Additional Welds | 15 |
| Cast Retention to Frame | 25 |
| Add Clasp to Frame (cast clasp) | 25 |
| Add Wrought Wire Clasp | 20 |
| Add Mesh Reinforcement | 25 |
| Add Tooth Colored Clasp | 30 |

| ADD'L REMOVABLE ITEN | MS |
|---|-----------|
| Bite Block with wax base plate (standard) | 10 |
| Bite Block with light-cured base plate (must note RX) | 15 |
| Bleaching Tray | 35 |
| Boxed Model | 5 |
| Boxed Model – Implant | 5 |
| CuSil – per gasket | 35 |
| Custom Tray – Light Cured | 15 |
| Digital Design - Denture (per arch) | 25 |
| Digital Design – Nightguard (per arch) | 10 |
| Duplicate Denture | 105 |
| Duplicate Model | 5 |
| Mounting | included |
| Patient ID (required in most states) | included |
| Printed Model – Full Arch (per arch) | 12 |
| Remove Teeth / Immediate Extraction (per tooth) | 5 |
| Reset | 5 |
| Teeth – Standard (per arch) | included |
| Teeth – Premium Ivoclar (per arch) | 58 |
| Teeth – Gold (Full / Open Face) approx per tooth | 150 - 250 |

Digital Design Turnaround: 24 to 48 hours

Turnaround Time: 8 business days | 6 business days (digital)

^{*}We believe in being open and honest with our Smile Partners. All the products noted on pages 2-4 are done by our offshore laboratory partner that we have used for over 15 years. For US made restorations we have an Elite option for crowns and our Prestige Removable Line (page 5).

PRESTIGE REMOVABLES (P)

Looking for a faster turnaround time? Try our Prestige Removables made in the US.

| FULL DENTURES | |
|--|-----|
| COMPLETE | 199 |
| includes mounting, set-up, Artic teeth, Lucitone 199 finish & patient ID | |
| Set-up (3 days) | 60 |
| Finish (4 days) | 79 |

| Ruch Available | for each stage. | add \$30 per dav |
|----------------|-----------------|------------------------------------|
| | | |

| METAL PARTIALS – excludes teeth | | |
|---|-----|--|
| COMPLETE | 259 | |
| includes frame, dup model, mounting, set-up, Lucitone 199 finish & patient ID | | |
| Add Flexible Clasp (per clasp) | 30 | |
| Frame only - SLM (selective laser melting) | 99 | |
| Set-up – Metal Partial (3 days) | 50 | |
| Finish – Metal Partial (4 days) | 75 | |
| SLM Frame: Digital - 5 Days Physical – 7 Da Rush service available: -1 day: \$115 / -2 Days: | | |

FLEXIBLE PARTIALS – excludes teeth

| includes duplicate model, mounting, set-up, & tcs® flexible finish | 149 |
|---|-----|
| Unilateral/Nesbit (2 teeth max) | 119 |
| Set-up – Flexible Partial (3 days) | 50 |
| Finish – Flexible Partial | 79 |
| Add Flexible clasp to frame (per clasp) | 30 |
| VisiClear Frame only (7 days / 9 days) | 149 |

Rush service available for Digital Cases: 3D Print Flex straight to finish: 4 Days: \$229 / 5 Days: \$205 VisiClear Frame only: 4 Days: \$199 / 5 Days: \$175

| ACRYLIC FLIPPERS – excludes teeth | |
|-----------------------------------|-----|
| Acrylic Flipper (1 to 4 teeth) | 100 |
| Acrylic Flipper (5+ teeth) | 115 |
| Add wrought wire clasps | 20 |
| 4 Days start to finish | |

4 Days start to finish

Rush service available (non-digital): add \$30 per day

| NIGHTGUARDS | |
|----------------------------|----|
| Soft (2 days) | 40 |
| Hard / Soft combo (3 days) | 50 |
| Hard – Acrylic (4 days) | 99 |

| ADD'L REMOVABLE | |
|---|----|
| Bite Block with light cure base plate (2 days) | 22 |
| Bleaching Tray (2 days) | 30 |
| Boxed Model | 10 |
| Custom Tray – Light Cured (2 days) | 22 |
| Digital Model – Full Arch (1 arch) | 12 |
| Duplicate Model | 10 |
| Mounting | 10 |
| Patient ID (required in most states) | 15 |
| Remove Teeth / Immediate Extraction (per tooth) | 5 |
| Reset (per arch, 1st is free) | 20 |
| Teeth – Standard per arch (Kulzer Artic) | 35 |
| Teeth – Premium per arch (Ivoclar BlueLine) | 86 |

| REPAIRS | |
|-------------------------|----|
| Acrylic Repair (2 days) | 59 |
| each additional tooth | 12 |
| Reline – Hard (2 days) | 70 |
| Cast Clasp (4 days) | 39 |
| Laser Weld (4 days) | 50 |
| Add Mesh Reinforcement | 20 |

TURNAROUND TIME

| Bite Block Custom Tray | 2 days |
|--------------------------------|-----------------|
| Setup | 4 days |
| Finish | 4 days |
| Flexible Finish | 6 – 8 days |
| Repairs | varies |
| SLM Frame (digital physical) | 5 days 7 days |

You must mark **PRESTIGE** on the Rx.

RETAINERS & ORTHDONTICS

| Acrylic Inlay 22 Acrylic Wrap 8 Adams Clasps 6 Anterior Bite Plane 11 Ball Clasps 5 Bands 6 C Clasps 7 Carve tooth off cast and prep site 11 Continuous Occlusal Reset 12 Custom Color or Design 15 Custom Palate 7 Decal added to acrylic palate 6 Diastema Spring 12 Essix Retainer (varies depending on complexity) 49 - 79 Finger Spring 6 Fixed Transpalatal (TPA) 45 Gelb Splint 117 Haas Expander 99 Hawley Retainer 49 Hawley Retainer w/ wrap around labial wire 87 Headgear Hooks 16 Indirect Bonding Trays 45 Invisible Retainer 40 Lingual 6x6 45 Nance/Rickinator 39 Palatal Piece 41 Pour Model - Ortho 13 | 8 working days in lab | |
|--|---|---------|
| Adams Clasps 6 Anterior Bite Plane 11 Ball Clasps 5 Bands 6 C Clasps 7 Carve tooth off cast and prep site 11 Continuous Occlusal Reset 12 Custom Color or Design 15 Custom Palate 7 Decal added to acrylic palate 6 Diastema Spring 12 Essix Retainer (varies depending on complexity) 49 - 79 Finger Spring 6 Fixed Transpalatal (TPA) 45 Gelb Splint 117 Haas Expander 99 Hawley Retainer 49 Hawley Retainer w/ wrap around labial wire 87 Headgear Hooks 16 Indirect Bonding Trays 45 Invisible Retainer 40 Lingual 6x6 45 Nance/Rickinator 39 Palatal Piece 41 Pontic Added 21 Pontic Added 21 Pour Model - Ortho 13 | Acrylic Inlay | 22 |
| Anterior Bite Plane 11 Ball Clasps 5 Bands 6 C Clasps 7 Carve tooth off cast and prep site 11 Continuous Occlusal Reset 12 Custom Color or Design 15 Custom Palate 7 Decal added to acrylic palate 6 Diastema Spring 12 Essix Retainer (varies depending on complexity) 49 - 79 Finger Spring 6 Fixed Transpalatal (TPA) 45 Gelb Splint 117 Haas Expander 99 Hawley Retainer 49 Hawley Retainer w/ wrap around labial wire 87 Headgear Hooks 16 Indirect Bonding Trays 45 Invisible Retainer 40 Lingual 6x6 45 Nance/Rickinator 39 Palatal Piece 41 Pontic Added 21 Pour Model - Ortho 13 Quad Helix 66 Rapid Palatal Expander (RPE) 70 <td>Acrylic Wrap</td> <td>8</td> | Acrylic Wrap | 8 |
| Ball Clasps 5 Bands 6 C Clasps 7 Carve tooth off cast and prep site 11 Continuous Occlusal Reset 12 Custom Color or Design 15 Custom Palate 7 Decal added to acrylic palate 6 Diastema Spring 12 Essix Retainer (varies depending on complexity) 49 - 79 Finger Spring 6 Fixed Transpalatal (TPA) 45 Gelb Splint 117 Haas Expander 99 Hawley Retainer 49 Hawley Retainer w/ wrap around labial wire 87 Headgear Hooks 16 Indirect Bonding Trays 45 Invisible Retainer 40 Lingual 6x6 45 Nance/Rickinator 39 Palatal Piece 41 Pontic Added 21 Pour Model - Ortho 13 Quad Helix 66 Rapid Palatal Expander (RPE) 70 Banded RPE 105 < | Adams Clasps | 6 |
| Bands 6 C Clasps 7 Carve tooth off cast and prep site 11 Continuous Occlusal Reset 12 Custom Color or Design 15 Custom Palate 7 Decal added to acrylic palate 6 Diastema Spring 12 Essix Retainer (varies depending on complexity) 49 - 79 Finger Spring 6 Fixed Transpalatal (TPA) 45 Gelb Splint 117 Haas Expander 99 Hawley Retainer 49 Hawley Retainer w/ wrap around labial wire 87 Headgear Hooks 16 Indirect Bonding Trays 45 Invisible Retainer 40 Lingual 6x6 45 Nance/Rickinator 39 Palatal Piece 41 Pontic Added 21 Pour Model - Ortho 13 Quad Helix 66 Rapid Palatal Expander (RPE) 70 Banded RPE 105 Reset per tooth 8 | Anterior Bite Plane | 11 |
| C Clasps 7 Carve tooth off cast and prep site 11 Continuous Occlusal Reset 12 Custom Color or Design 15 Custom Palate 7 Decal added to acrylic palate 6 Diastema Spring 12 Essix Retainer (varies depending on complexity) 49 - 79 Finger Spring 6 Fixed Transpalatal (TPA) 45 Gelb Splint 117 Haas Expander 99 Hawley Retainer 49 Hawley Retainer w/ wrap around labial wire 87 Headgear Hooks 16 Indirect Bonding Trays 45 Invisible Retainer 40 Lingual 6x6 45 Nance/Rickinator 39 Palatal Piece 41 Pontic Added 21 Pour Model - Ortho 13 Quad Helix 66 Rapid Palatal Expander (RPE) 70 Banded RPE 105 Reset per tooth 8 Rickets Labial Wire 17 </td <td>Ball Clasps</td> <td>5</td> | Ball Clasps | 5 |
| Carve tooth off cast and prep site 11 Continuous Occlusal Reset 12 Custom Color or Design 15 Custom Palate 7 Decal added to acrylic palate 6 Diastema Spring 12 Essix Retainer (varies depending on complexity) 49 - 79 Finger Spring 6 Fixed Transpalatal (TPA) 45 Gelb Splint 117 Haas Expander 99 Hawley Retainer 49 Hawley Retainer w/ wrap around labial wire 87 Headgear Hooks 16 Indirect Bonding Trays 45 Invisible Retainer 40 Lingual 6x6 45 Nance/Rickinator 39 Palatal Piece 41 Pontic Added 21 Pour Model - Ortho 13 Quad Helix 66 Rapid Palatal Expander (RPE) 70 Banded RPE 105 Reset per tooth 8 Rickets Labial Wire 17 Rollo Bands (each) | Bands | 6 |
| Continuous Occlusal Reset 12 Custom Color or Design 15 Custom Palate 7 Decal added to acrylic palate 6 Diastema Spring 12 Essix Retainer (varies depending on complexity) 49 - 79 Finger Spring 6 Fixed Transpalatal (TPA) 45 Gelb Splint 117 Haas Expander 99 Hawley Retainer 49 Hawley Retainer w/ wrap around labial wire 87 Headgear Hooks 16 Indirect Bonding Trays 45 Invisible Retainer 40 Lingual 6x6 45 Nance/Rickinator 39 Palatal Piece 41 Pontic Added 21 Pour Model - Ortho 13 Quad Helix 66 Rapid Palatal Expander (RPE) 70 Banded RPE 105 Reset per tooth 8 Rickets Labial Wire 17 Rollo Bands (each) 11 Soldered Spring 25 | C Clasps | 7 |
| Custom Color or Design 15 Custom Palate 7 Decal added to acrylic palate 6 Diastema Spring 12 Essix Retainer (varies depending on complexity) 49 - 79 Finger Spring 6 Fixed Transpalatal (TPA) 45 Gelb Splint 117 Haas Expander 99 Hawley Retainer 49 Hawley Retainer w/ wrap around labial wire 87 Headgear Hooks 16 Indirect Bonding Trays 45 Invisible Retainer 40 Lingual 6x6 45 Nance/Rickinator 39 Palatal Piece 41 Pontic Added 21 Pour Model - Ortho 13 Quad Helix 66 Rapid Palatal Expander (RPE) 70 Banded RPE 105 Reset per tooth 8 Rickets Labial Wire 17 Rollo Bands (each) 11 Soldered Spring 25 Space Maintainer 29 - 49 <td>Carve tooth off cast and prep site</td> <td>11</td> | Carve tooth off cast and prep site | 11 |
| Custom Palate 7 Decal added to acrylic palate 6 Diastema Spring 12 Essix Retainer (varies depending on complexity) 49 - 79 Finger Spring 6 Fixed Transpalatal (TPA) 45 Gelb Splint 117 Haas Expander 99 Hawley Retainer 49 Hawley Retainer w/ wrap around labial wire 87 Headgear Hooks 16 Indirect Bonding Trays 45 Invisible Retainer 40 Lingual 6x6 45 Nance/Rickinator 39 Palatal Piece 41 Pontic Added 21 Pour Model - Ortho 13 Quad Helix 66 Rapid Palatal Expander (RPE) 70 Banded RPE 105 Reset per tooth 8 Rickets Labial Wire 17 Rollo Bands (each) 11 Soldered Clasps 14 Soldered Spring 25 Space Maintainer 29 - 49 < | Continuous Occlusal Reset | 12 |
| Decal added to acrylic palate 6 Diastema Spring 12 Essix Retainer (varies depending on complexity) 49 - 79 Finger Spring 6 Fixed Transpalatal (TPA) 45 Gelb Splint 117 Haas Expander 99 Hawley Retainer 49 Hawley Retainer w/ wrap around labial wire 87 Headgear Hooks 16 Indirect Bonding Trays 45 Invisible Retainer 40 Lingual 6x6 45 Nance/Rickinator 39 Palatal Piece 41 Pontic Added 21 Pour Model - Ortho 13 Quad Helix 66 Rapid Palatal Expander (RPE) 70 Banded RPE 105 Reset per tooth 8 Rickets Labial Wire 17 Rollo Bands (each) 11 Soldered Clasps 14 Soldered Spring 25 Space Maintainer 29 - 49 Support Wires 5 <td>Custom Color or Design</td> <td>15</td> | Custom Color or Design | 15 |
| Diastema Spring 12 Essix Retainer (varies depending on complexity) 49 - 79 Finger Spring 6 Fixed Transpalatal (TPA) 45 Gelb Splint 117 Haas Expander 99 Hawley Retainer 49 Hawley Retainer w/ wrap around labial wire 87 Headgear Hooks 16 Indirect Bonding Trays 45 Invisible Retainer 40 Lingual 6x6 45 Nance/Rickinator 39 Palatal Piece 41 Pontic Added 21 Pour Model - Ortho 13 Quad Helix 66 Rapid Palatal Expander (RPE) 70 Banded RPE 105 Reset per tooth 8 Rickets Labial Wire 17 Rollo Bands (each) 11 Soldered Clasps 14 Soldered Spring 25 Space Maintainer 29 - 49 Support Wires 5 | Custom Palate | 7 |
| Essix Retainer (varies depending on complexity) 49 - 79 Finger Spring 6 Fixed Transpalatal (TPA) 45 Gelb Splint 117 Haas Expander 99 Hawley Retainer 49 Hawley Retainer w/ wrap around labial wire 87 Headgear Hooks 16 Indirect Bonding Trays 45 Invisible Retainer 40 Lingual 6x6 45 Nance/Rickinator 39 Palatal Piece 41 Pontic Added 21 Pour Model - Ortho 13 Quad Helix 66 Rapid Palatal Expander (RPE) 70 Banded RPE 105 Reset per tooth 8 Rickets Labial Wire 17 Rollo Bands (each) 11 Soldered Clasps 14 Soldered Spring 25 Space Maintainer 29 - 49 Support Wires 5 | Decal added to acrylic palate | 6 |
| Finger Spring 6 Fixed Transpalatal (TPA) 45 Gelb Splint 117 Haas Expander 99 Hawley Retainer 49 Hawley Retainer w/ wrap around labial wire 87 Headgear Hooks 16 Indirect Bonding Trays 45 Invisible Retainer 40 Lingual 6x6 45 Nance/Rickinator 39 Palatal Piece 41 Pontic Added 21 Pour Model - Ortho 13 Quad Helix 66 Rapid Palatal Expander (RPE) 70 Banded RPE 105 Reset per tooth 8 Rickets Labial Wire 17 Rollo Bands (each) 11 Soldered Clasps 14 Soldered Spring 25 Space Maintainer 29 - 49 Support Wires 5 | Diastema Spring | 12 |
| Fixed Transpalatal (TPA) 45 Gelb Splint 117 Haas Expander 99 Hawley Retainer 49 Hawley Retainer w/ wrap around labial wire 87 Headgear Hooks 16 Indirect Bonding Trays 45 Invisible Retainer 40 Lingual 6x6 45 Nance/Rickinator 39 Palatal Piece 41 Pontic Added 21 Pour Model - Ortho 13 Quad Helix 66 Rapid Palatal Expander (RPE) 70 Banded RPE 105 Reset per tooth 8 Rickets Labial Wire 17 Rollo Bands (each) 11 Soldered Clasps 14 Soldered Spring 25 Space Maintainer 29 - 49 Support Wires 5 | Essix Retainer (varies depending on complexity) | 49 - 79 |
| Gelb Splint 117 Haas Expander 99 Hawley Retainer 49 Hawley Retainer w/ wrap around labial wire 87 Headgear Hooks 16 Indirect Bonding Trays 45 Invisible Retainer 40 Lingual 6x6 45 Nance/Rickinator 39 Palatal Piece 41 Pontic Added 21 Pour Model - Ortho 13 Quad Helix 66 Rapid Palatal Expander (RPE) 70 Banded RPE 105 Reset per tooth 8 Rickets Labial Wire 17 Rollo Bands (each) 11 Soldered Clasps 14 Soldered Spring 25 Space Maintainer 29 - 49 Support Wires 5 | Finger Spring | 6 |
| Haas Expander99Hawley Retainer49Hawley Retainer w/ wrap around labial wire87Headgear Hooks16Indirect Bonding Trays45Invisible Retainer40Lingual 6x645Nance/Rickinator39Palatal Piece41Pontic Added21Pour Model - Ortho13Quad Helix66Rapid Palatal Expander (RPE)70Banded RPE105Reset per tooth8Rickets Labial Wire17Rollo Bands (each)11Soldered Clasps14Soldered Spring25Space Maintainer29 - 49Support Wires5 | Fixed Transpalatal (TPA) | 45 |
| Hawley Retainer49Hawley Retainer w/ wrap around labial wire87Headgear Hooks16Indirect Bonding Trays45Invisible Retainer40Lingual 6x645Nance/Rickinator39Palatal Piece41Pontic Added21Pour Model - Ortho13Quad Helix66Rapid Palatal Expander (RPE)70Banded RPE105Reset per tooth8Rickets Labial Wire17Rollo Bands (each)11Soldered Clasps14Soldered Spring25Space Maintainer29 - 49Support Wires5 | Gelb Splint | 117 |
| Hawley Retainer w/ wrap around labial wire Headgear Hooks 16 Indirect Bonding Trays 45 Invisible Retainer 40 Lingual 6x6 Nance/Rickinator Palatal Piece 41 Pontic Added 21 Pour Model - Ortho 13 Quad Helix 66 Rapid Palatal Expander (RPE) Reset per tooth Rickets Labial Wire 70 Rollo Bands (each) Soldered Clasps Space Maintainer 29 - 49 Support Wires 5 | Haas Expander | 99 |
| Headgear Hooks 16 Indirect Bonding Trays 45 Invisible Retainer 40 Lingual 6x6 45 Nance/Rickinator 39 Palatal Piece 41 Pontic Added 21 Pour Model - Ortho 13 Quad Helix 66 Rapid Palatal Expander (RPE) 70 Banded RPE 105 Reset per tooth 8 Rickets Labial Wire 17 Rollo Bands (each) 11 Soldered Clasps 14 Soldered Spring 25 Space Maintainer 29 - 49 Support Wires 5 | Hawley Retainer | 49 |
| Indirect Bonding Trays 45 Invisible Retainer 40 Lingual 6x6 45 Nance/Rickinator 39 Palatal Piece 41 Pontic Added 21 Pour Model - Ortho 13 Quad Helix 66 Rapid Palatal Expander (RPE) 70 Banded RPE 105 Reset per tooth 8 Rickets Labial Wire 17 Rollo Bands (each) 11 Soldered Clasps 14 Soldered Spring 25 Space Maintainer 29 - 49 Support Wires 5 | Hawley Retainer w/ wrap around labial wire | 87 |
| Invisible Retainer 40 Lingual 6x6 45 Nance/Rickinator 39 Palatal Piece 41 Pontic Added 21 Pour Model - Ortho 13 Quad Helix 66 Rapid Palatal Expander (RPE) 70 Banded RPE 105 Reset per tooth 8 Rickets Labial Wire 17 Rollo Bands (each) 11 Soldered Clasps 14 Soldered Spring 25 Space Maintainer 29 - 49 Support Wires 5 | Headgear Hooks | 16 |
| Lingual 6x6 45 Nance/Rickinator 39 Palatal Piece 41 Pontic Added 21 Pour Model - Ortho 13 Quad Helix 66 Rapid Palatal Expander (RPE) 70 Banded RPE 105 Reset per tooth 8 Rickets Labial Wire 17 Rollo Bands (each) 11 Soldered Clasps 14 Soldered Spring 25 Space Maintainer 29 - 49 Support Wires 5 | Indirect Bonding Trays | 45 |
| Nance/Rickinator 39 Palatal Piece 41 Pontic Added 21 Pour Model - Ortho 13 Quad Helix 66 Rapid Palatal Expander (RPE) 70 Banded RPE 105 Reset per tooth 8 Rickets Labial Wire 17 Rollo Bands (each) 11 Soldered Clasps 14 Soldered Spring 25 Space Maintainer 29 - 49 Support Wires 5 | Invisible Retainer | 40 |
| Palatal Piece 41 Pontic Added 21 Pour Model - Ortho 13 Quad Helix 66 Rapid Palatal Expander (RPE) 70 Banded RPE 105 Reset per tooth 8 Rickets Labial Wire 17 Rollo Bands (each) 11 Soldered Clasps 14 Soldered Spring 25 Space Maintainer 29 - 49 Support Wires 5 | Lingual 6x6 | 45 |
| Pontic Added 21 Pour Model - Ortho 13 Quad Helix 66 Rapid Palatal Expander (RPE) 70 Banded RPE 105 Reset per tooth 8 Rickets Labial Wire 17 Rollo Bands (each) 11 Soldered Clasps 14 Soldered Spring 25 Space Maintainer 29 - 49 Support Wires 5 | Nance/Rickinator | 39 |
| Pour Model - Ortho 13 Quad Helix 66 Rapid Palatal Expander (RPE) 70 Banded RPE 105 Reset per tooth 8 Rickets Labial Wire 17 Rollo Bands (each) 11 Soldered Clasps 14 Soldered Spring 25 Space Maintainer 29 - 49 Support Wires 5 | Palatal Piece | 41 |
| Quad Helix 66 Rapid Palatal Expander (RPE) 70 Banded RPE 105 Reset per tooth 8 Rickets Labial Wire 17 Rollo Bands (each) 11 Soldered Clasps 14 Soldered Spring 25 Space Maintainer 29 - 49 Support Wires 5 | Pontic Added | 21 |
| Rapid Palatal Expander (RPE) 70 Banded RPE 105 Reset per tooth 8 Rickets Labial Wire 17 Rollo Bands (each) 11 Soldered Clasps 14 Soldered Spring 25 Space Maintainer 29 - 49 Support Wires 5 | Pour Model - Ortho | 13 |
| Banded RPE 105 Reset per tooth 8 Rickets Labial Wire 17 Rollo Bands (each) 11 Soldered Clasps 14 Soldered Spring 25 Space Maintainer 29 - 49 Support Wires 5 | Quad Helix | 66 |
| Reset per tooth 8 Rickets Labial Wire 17 Rollo Bands (each) 11 Soldered Clasps 14 Soldered Spring 25 Space Maintainer 29 - 49 Support Wires 5 | Rapid Palatal Expander (RPE) | 70 |
| Rickets Labial Wire 17 Rollo Bands (each) 11 Soldered Clasps 14 Soldered Spring 25 Space Maintainer 29 - 49 Support Wires 5 | Banded RPE | 105 |
| Rollo Bands (each)11Soldered Clasps14Soldered Spring25Space Maintainer29 - 49Support Wires5 | Reset per tooth | 8 |
| Soldered Clasps14Soldered Spring25Space Maintainer29 - 49Support Wires5 | Rickets Labial Wire | 17 |
| Soldered Spring25Space Maintainer29 - 49Support Wires5 | Rollo Bands (each) | 11 |
| Space Maintainer29 - 49Support Wires5 | Soldered Clasps | 14 |
| Support Wires 5 | Soldered Spring | 25 |
| | Space Maintainer | 29 - 49 |
| Thumb/Tongue Crib 75 | Support Wires | 5 |
| | Thumb/Tongue Crib | 75 |

SNORING & SLEEP APNEA*

| 7 working days in lab | |
|---|-----|
| Silent Nite® Sleep Appliance* | 215 |
| OASYS Hinge Appliance™ sleep apnea device | 459 |
| EMA® sleep apnea device | 289 |
| TAP 3 TL sleep apnea device | 459 |
| dreamTAP® sleep apnea device | 559 |

^{*}Buy 2 appliances and take \$25 off the second appliance

SMILE SHAPERS – CLEAR ALIGNERS

Fabricated out of Zendura, a biocompatible thermoplastic, Smile Shapers are virtually invisible, durable and comfortable. All aligners are expertly designed and created to ensure their precision and durability.

| Clear Aligner (per aligner) - 7 days | 35 |
|---|-----|
| Treatment Plan / Clinical Proposal – 3 days | 85 |
| Scan Fee (when physical models sent) | 30 |
| Attachment Template | 25 |
| Extended Care Package | 195 |
| Replacement Aligner | 30 |
| Final Retainer | 35 |
| Starter Kit | 99 |

SPORTS MOUTHGUARDS*

Our mouthguards can be customized for youth, amateur and even professional athletes and teams. Call us for bulk ordering and quotes including adding your team's colors and logos.

| Basic Sports Mouthguard | 75 |
|---|----|
| Add Helmet Strap | 15 |
| Add Name | 5 |
| *Pricing will vary based on thickness of appliance. | |

| RESTORATION / CASE TYPE | WORKING DAYS |
|-------------------------|-----------------|
| Crown & Bridge | 8 |
| from digital file | 6 |
| Removables (each stage) | 8 |
| from digital file | 6 |
| Prestige Removables | Varies |
| Orthodontics | 8 |
| Digital Design | 24 – 48 hours |

Please allow full working time for each product selected.

Working times are <u>NOT</u> guaranteed and do <u>NOT</u> include weekends or holidays or shipping times.

Time of pickup and delivery may affect turnaround time.

Rush service available on select restorations but must be prescheduled. To inquire or preschedule rush service, please call 877.622.3533 before the case is shipped.

Hours of Operation

Monday - Friday: 8:00 a.m. to 4:30 p.m. EST

Digital File Cut-Off Time: 4:00 PM (EST)

Prices shown are per unit, per arch or per appliance and do not include applicable taxes.

Prices and product availability subject to change without notice.

Volume discounts available, inquire within.

Simply Your Experience with DDX

With our online portal you can save valuable time and eliminate unneeded phone calls. With DDX, you can:

- Schedule case pickups
- View your case's active status
- Track your shipped cases
- Pay your lab bill
- Setup auto payments
- View price lists
- Send digital impression files and patient photos
- View your balance and statements

Send Us Your Next Case

- Fill out a Prescription Rx You may download any Rx at mabeldental.com/downloads/ or by calling 877.622.3533
- Call 877.622.3533 or email: support@mabeldental.com to establish your account
- Shipping Supplies Call 877.622.3533 or send an email to: support@mabeldental.com to request shipping supplies including boxes as needed. You can print UPS labels directly from our website: www.mabeldental.com/print-ups-label/
- 4. Package and Ship Your Case. Separate impressions, models and Rx into bags to ensure accurate identification. Place protective foam or cushion over cases, seal box securely with packing tape, and attach mailing label. Save money by combining multiple cases in one box.

UPS Shipping: Give your case to any UPS Driver, UPS Drop-Box or UPS store. Or call us to schedule a UPS pickup* 877.622.3533, a small fee applies.

Courier: Call 877.622.3533 to schedule an in-office case pickup with our courier in most areas within Ohio

 Digital Cases: upload directly via your DDX account, email files to support@mabeldental.com or upload via our File Upload webpage

Shipping: \$9* / each way

*pricing varies depending on location and timing

If you ship to us using your shipping label, we will ship back to you at no charge. Digital cases will incur an outbound shipping charge.

Terms: All accounts must be signed up for credit card auto payment or include a check with each case.

Accounts not paid within the stated terms will be subject to COD status and a late charge of 2 percent of the unpaid balance. Cases held until past due balance is paid in full.

Limited Warranty: For warranty terms and conditions of limitation of liability, visit our website: mabeldental.com/policies-and-warranty